
IFU 03 01 0 00 REV 4_02-2026

3DIEMME S.r.l.
Via Risorgimento 9
22063 Cantù (CO) - Italia
www.3diemme.it
support@3diemme.it
+39 0317073353

Contents

(EN) English - INSTRUCTIONS FOR USE	2
(IT) Italiano - ISTRUZIONI D'USO	5
(ES) Español - INSTRUCCIONES DE USO	8
(FR) Français - MODE D'EMPLOI	11
(BG) Български - ИНСТРУКЦИИ ЗА УПОТРЕБА	17
(CS) Čeština - NÁVOD K POUŽITÍ	20
(DA) Dansk - BRUGSANVISNING	23
(DE) Deutsch - GEBRAUCHSANLEITUNG	25
(HU) Magyar - HASZNÁLATI UTASÍTÁS	28
(PT) INSTRUÇÕES PARA O USO	31
(RO) Română - INSTRUȚIUNI DE UTILIZARE	34
(SK) Slovenčina - NÁVOD NA POUŽITIE	37
(SL) Slovenščina - NAVODILA ZA UPORABO	40
(SV) Suomi - BRUKSANVISNING	43

If any serious incident and/ or complaint occurs using a sleeve, please contact 3DIEMME or your distributor accordingly. Our company will conduct a proper analysis procedure and will report the event when applicable to the competent authority of the Member State in which the user and/or patient is established.

Trademark Notice

In certain jurisdictions, including the United States and the European Union, RealGUIDE is a registered trademark of 3Diemme S.r.l.. All references of the symbol ® are limited to the relevant jurisdictions where the trademark is registered.

1. Product identification

RealGUIDE Sleeve is a metal ring that is manually inserted into a 3D printed surgical guide, specifically designed for a Patient and based on the Doctor's project made in RealGUIDE software.

2. Materials

RealGUIDE Sleeves are made of stainless steel KLEINOX 4598 or in Titan grade 5 TiAl6V4 ELI or PEEK.

3. Intended Use

The *RealGUIDE Instruments and Accessories* are intended to contribute in the creation of the implant site, the insertion and extraction of the dental implant and for the connection and disconnection of prosthetic components to the dental implant.

4. Indication for use

RealGUIDE Sleeves are intended to:

- Drive in depth and inclination the surgical drills to prepare the bone for the dental implant insertion;
- Drive in depth, inclination and phase the implant mounts, connected to the implants through a fixation screw, so that the dental implants, virtually placed in the software, are transferred correctly into the Patient's mouth;
- Drive in depth and inclination the eventual fixation tools used by the surgeon to keep the surgical guide in position in case of totally edentulous patients.

5. Contraindications

The use of RealGUIDE Instruments and Accessories has no contraindications if used for the purpose intended by the manufacturer.

The use of dental guided surgery instruments is contra-indicated where dental guided surgery and thus implant placement is contra-indicated, such as in vulnerable patients (e.g. pregnant women, pediatric and elderly patients).

6. Intended users

The handling and use of the RealGUIDE Instruments and Accessories must be carried out by qualified professional personnel (implantologist, dentist or maxillofacial surgeon) who has the necessary qualifications required according

to the legislation in force in each country. It is necessary to have adequate and specialised training on the guided surgery technique. The intervention must be performed in a protected environment which is considered adequate from a hygienic point of view and with appropriate instrumentation.

7. Target patient population

The device is intended for all patients requiring dental guided surgery, both males or females of any age.

8. Storage

The goods are sent in a protective packaging in order to prevent damage being caused during storage. Store the products at room temperature.

9. Cleaning and disinfection Instructions

Once inserted in the surgical guide, sleeves are cleaned together with the guide.

10. Methods for Use

Glue the sleeve (**check the sleeve compatibility with the planned implant before performing surgery**) in the surgical guide hole paying attention to the profile orientation (if available).

11. Warnings

- RealGUIDE Sleeves are supplied non-sterile and must be cleaned before use (see Cleaning Instructions). Using a surgical instrument not adequately cleaned exposes doctors, operators and patients to serious risk of infection and cross-infection.
- This product is for single use and doesn't need maintenance. In fact, the use of the surgical guide could deform the sleeve's hole (if drills are not inserted on axis with the sleeve).

12. Disposal

After terminating the procedure, or in case the device needs to be disposed of, it must be treated as biological waste and disposal must follow the clinic's requirements and disposal instructions.


13. Symbol glossary

	Caution		Prescription only (USA market)
 eIFU	Consult instructions for use		Batch Code
	CE Marking		Manufacturer
	Non-sterile		Medical Device
	Catalogue Number		Do not reuse

	Do not use if package is damaged		Swiss Authorised Representative
---	-------------------------------------	---	------------------------------------

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